

Best.-Nr.: 2030019066  
FAID0005

DE  
EN  
PL

**EN**

3

.....  
Please refer to the graphics in the German Installation and Operating Instructions.

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Emergency shower with pull rod

Please refer to the graphics in the German installation and operating instructions.

## Table of contents

1. Abbreviations and units . . . . .	3
2. Key . . . . .	4
3. Warranty . . . . .	4
4. Important notes . . . . .	4
<b>Description of product</b>	
5. Application . . . . .	4
6. Technical specifications . . . . .	4
7. Scope of delivery . . . . .	5
8. Dimensions . . . . .	5
9. Installation example . . . . .	5
<b>Assembly, function and commissioning</b>	
10. Assembly for actuation from right . . . . .	5
11. Assembly for actuation from left . . . . .	7
12. Function . . . . .	8
<b>Maintenance</b>	
13. Maintenance and care . . . . .	8
14. Troubleshooting . . . . .	8
15. Spare parts . . . . .	8

## 1. Abbreviations and units

Best.-Nr.                                      Aquarotter order number

Conversion                                      1 mm = 0.03937 inches

1 inch = 25.4 mm

All length specifications in the graphics are in mm.

## 2. Key

---

### Warning!

Failure to observe can result in bodily injury or even death.

### Caution!

Failure to observe can result in material damage.

### Important!

Failure to observe can cause the product to malfunction.

Useful information for optimum handling of the product.

## 3. Warranty

---

Liability is accepted in accordance with the General Terms and Conditions of Business and Supply.

Use only genuine spare parts!

## 4. Important notes

---

- Installation, commissioning and maintenance must be performed only by a qualified expert in accordance with the instructions provided, legal requirements and recognised engineering standards.
- All technical connection regulations specified by the local water and electricity supply companies must be observed.
- We reserve the right to make changes.

## 5. Application

---

Emergency showers are prescribed first aid installations for work stations where hazardous materials are being handled.

## 6. Technical specifications

---

Minimum flow pressure: 0.5 bar

Maximum stagnation pressure: 10 bar

Volumetric flow:

0.5 bar flow pressure: 45 l/m

1.0 bar flow pressure: 65 l/m

1.5 bar flow pressure: 80 l/m

2.0 bar flow pressure: 90 l/m

2.5 bar flow pressure: 100 l/m

3.0 bar flow pressure: 110 l/m

Assembly height: Lower edge of shower head (220 ± 10) cm above floor as per EN 15154 Part 1

## 7. Scope of delivery

---

Pos.	Pcs.	Name
1	1	Shower head
2	1	Wall bracket
3	1	Wall panel
4	1	Ball valve
5	1	Wall support
6	1	Pipe extension*
7	1	Connecting pipe
8	1	Hand actuator with plugging option
9	1	Accessories kit: First aid sign
	1	Installation and operating instructions

\* Not needed for actuation on the left.

## 8. Dimensions

---

A: Design with pull rod on right

B: Design with pull rod on left

C: For special designs, see the customer drawing for dimensions

## 9. Installation example

---

## 10. Assembly for actuation from right

---

### Important!

- Prior to the assembly process, flush the pipework in accordance with DIN 1988.
- Mount the emergency shower in accordance with EN 15154 Part 1.

 The manual actuator of the emergency shower can be plugged.

 To protect the tiles, use masking tape when marking and drilling. Drill at low speed.

 Depending on the composition of the wall, special wall plugs may have to be used (to be provided by the customer).

**10.1** Screw the connecting pipe (b) into the wall plate (a) so that a seal is formed.

**10.2** Screw the pipe extension (c) onto the connecting pipe (b) with an Allen key to form a seal.

- 10.3** Screw in the closed ball valve (d) as far as possible, creating a seal and bringing it into position as follows:
- The valve axis must point forwards (and not to the wall).
  - The pivot motion of the valve lever must be aligned so that it is parallel to the wall.
- 10.4** Screw the wall support (e) onto the ball valve (d) to form a seal and align flush with the wall plate (a).
-  For a door leaf width from 1 m, check whether the door frame obstructs the assembly of the pull rod or wall support.
- 10.5** Door leaf width up to 1 m:  
Align the wall panel (a) centrally above the door,  $2310 \pm 100$  mm above the finished floor in accordance with the water intake, and mark the drill holes.
- 10.6** Align the wall support (e) in accordance with the dimensions and mark the drill holes.
- 10.7** Drill in accordance with these markings.
- 10.8** Insert the wall plugs.
- 10.9** Screw on the wall plate (a) and the wall support (e).
- 10.10** Place the wall support (f) onto the guide (h) on the manual actuator and secure with the grub screw (g).
- 10.11** Place the valve lever into the "close" position.  
A: Closed  
B: Open
- 10.12** Screw the manual actuator (k) into the fork head (j) of the deflection lever and lock with a nut.
- 10.13** Using the spirit level, align the manual actuator (k) so that it is vertical.
- 10.14** Mark drill holes for the wall support (l) 40 mm above the triangular handle.
- 10.15** Drill in accordance with these markings.
- 10.16** Insert the wall plugs.
- 10.17** Screw on the wall support.
- 10.18** Mount the wall arm in a sealing manner and bring it into position.
- 10.19** Screw the shower head (n) onto the wall arm (m) so that it seals.
- 10.20** Establish the water connection.
- 10.21** Check the water connection for leaks.
- 10.22** Plug the manual actuator if need be.

## 11. Assembly for actuation from left

---

### Important!

- Prior to the assembly process, flush the pipework in accordance with DIN 1988.
- Mount the emergency shower in accordance with EN 15154 Part 1.

 The manual actuator of the emergency shower can be plugged.

 To protect the tiles, use masking tape when marking and drilling. Drill at low speed.

 Depending on the composition of the wall, special wall plugs may have to be used (to be provided by the customer).

**11.1** Screw the connecting pipe (a) into the wall plate (b) so that a seal is formed.

**11.2** Screw in the closed ball valve (c) as far as possible, creating a seal and bringing it into position as follows:

- The valve axis must point forwards (and not to the wall).
- The pivot motion of the valve lever must be aligned so that it is parallel to the wall.

**11.3** Screw the wall support (d) onto the ball valve (c) to form a seal and align flush with the wall plate (b).

 For a door leaf width from 1 m, check whether the door frame obstructs the assembly of the pull rod or wall support.

**11.4** Door leaf width up to 1 m:

Align the wall panel (b) centrally above the door,  $2310 \pm 100$  mm above the finished floor in accordance with the water intake, and mark the drill holes.

**11.5** Align the wall support (d) in accordance with the dimensions and mark the drill holes.

**11.6** Drill in accordance with these markings.

**11.7** Insert the wall plugs.

**11.8** Screw on the wall plate (b) and the wall support (d).

**11.9** Place the wall support (f) onto the guide (h) on the manual actuator and secure with the grub screw (g).

**11.10** Place the valve lever into the "close" position.

A: Closed

B: Open

**11.11** Screw the manual actuator (j) into the fork head (h) of the deflection lever and lock with a nut.

**11.12** Using the spirit level, align the manual actuator (j) so that it is vertical.

**11.13** Mark drill holes for the wall support (k) 40 mm above the triangular handle.

**11.14** Drill in accordance with these markings.

**11.15** Insert the wall plugs.

**11.16** Screw on the wall support.

**11.17** Screw the wall bracket (l) into the shower head (m) so that a seal is formed.

**11.18** Seal in the shower head with the wall bracket and bring into position.

- 11.19 Establish the water connection.
- 11.20 Check the water connection for leaks.
- 11.21 Plug the manual actuator if need be.

## 12. Function

---

12.1 Pull the manual actuator.

- The water will flow.

 The valve is not self-closing!

12.2 Place the manual actuator back in its start position.

- The water stops flowing.

## 13. Maintenance and care

---

Check the emergency shower once a month to make sure that it is in working order. (Guidelines for Laboratories BGI/GUV-I 850-0)

As a measure for the prevention of microbiological contamination, it is recommended that you replace the fill water in the fitting at shorter intervals by opening the valve.

Use suitable cleaning agents which do not attack the tap in a proper manner and rinse with water after use. High-pressure cleaners must not be used for cleaning.

## 14. Troubleshooting

---

Fault	Cause	Remedy
Water does not flow	– Water supply interrupted	⇒ Restore it
Manual actuator is catching	– Ball valve not correctly aligned	⇒ Align it
	– Pull rod or manual actuator bent	⇒ Straighten it
	– Ball valve stiff	⇒ Replace it
Body shower is dripping	– Ball valve defective	⇒ Replace it

If you cannot rectify a fault or if it is not in the list, please contact our customer service department!

## 15. Spare parts

---

Description	Order No.
Ball valve .....	2030025020

 Other spare parts on request; see the wall bracket for material numbers.

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KWC ME LLC Ras Al Khaimah,
KWC DVS Ltd - Northern Office
KWC DVS Ltd - Southern Office
KWC Aquarotter GmbH
KWC Nordics Oy
**KWC Austria GmbH**